

PF-04628935

Chemical Properties


CAS No. : 1383719-97-8

Formula: C₂₄H₂₆ClN₇O₅

Molecular Weight: 496.03

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PF-04628935 is a potent antagonist/inverse agonist of the ghrelin receptor, growth hormone secretagogue receptor 1a (GHS-R1a) with an IC ₅₀ of 4.6 nM. PF-04628935 is orally bioavailable and brain penetrant.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.016 mL	10.080 mL	20.1601 mL
5 mM	0.4032 mL	2.016 mL	4.032 mL
10 mM	0.2016 mL	1.008 mL	2.016 mL
50 mM	0.0403 mL	0.2016 mL	0.4032 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481