


RH01617

Chemical Properties

CAS No. : 302333-18-2
Formula: C₂₂H₂₆N₂O₅S
Molecular Weight: 430.52
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	RH01617 possesses potent inhibitory activities against Kv1.5.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3228 mL	11.6139 mL	23.2277 mL
5 mM	0.4646 mL	2.3228 mL	4.6455 mL
10 mM	0.2323 mL	1.1614 mL	2.3228 mL
50 mM	0.0465 mL	0.2323 mL	0.4646 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Guo X, Yang Q, Xu J, Zhang L, Chu H, Yu P, Zhu Y, Wei J, Chen W, Zhang Y, Zhang X, Sun H, Tang Y, You Q. Design and bio-evaluation of indole derivatives as potent Kv1.5 inhibitors. Bioorg Med Chem. 2013 Nov 1;21(21):6466-76.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481