


Ro 31-2201

Chemical Properties

CAS No. : 88851-65-4
Formula: C₂₀H₂₅N₃O₆
Molecular Weight: 403.43
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Ro 31-2201 is a potent inhibitor of angiotensin converting enzyme.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4787 mL	12.3937 mL	24.7874 mL
5 mM	0.4957 mL	2.4787 mL	4.9575 mL
10 mM	0.2479 mL	1.2394 mL	2.4787 mL
50 mM	0.0496 mL	0.2479 mL	0.4957 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Attwood MR, Francis RJ, Hassall CH, Kröhn A, Lawton G, Natoff IL, Nixon JS, Redshaw S, Thomas WA. New potent inhibitors of angiotensin converting enzyme. FEBS Lett. 1984 Jan 9;165(2):201-6. PubMed PMID: 6319181.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481