


Saviprazole

Chemical Properties

CAS No. : 121617-11-6
 Formula: C₁₅H₁₀F₇N₃O₂S₂
 Molecular Weight: 461.38
 Appearance: no data available
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Saviprazole (Hoe-731), a proton pump inhibitor, is used potentially for treatment of gastric ulcer.
Targets(IC50)	Proton pump

Solubility Information

Solubility	DMSO: 5 mg/ml (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1674 mL	10.8371 mL	21.6741 mL
5 mM	0.4335 mL	2.1674 mL	4.3348 mL
10 mM	0.2167 mL	1.0837 mL	2.1674 mL
50 mM	0.0433 mL	0.2167 mL	0.4335 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Herling AW, et al. Gastric acid inhibitory profile of saviprazole (HOE 731) compared to omeprazole. Pharmacology. 1991;43(6):293-303.

Inhibitor • Natural Compounds • Compound Libraries • Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481