


SCH-23982

Chemical Properties

CAS No. : 99234-87-4
Formula: C₁₇H₁₉NO
Molecular Weight: 253.34
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	SCH-23982 is a dopamine D1 receptor antagonist.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9473 mL	19.7363 mL	39.4726 mL
5 mM	0.7895 mL	3.9473 mL	7.8945 mL
10 mM	0.3947 mL	1.9736 mL	3.9473 mL
50 mM	0.0789 mL	0.3947 mL	0.7895 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Lillrank SM, Lipska BK, Weinberger DR, Fredholm BB, Fuxe K, Ferré S. Adenosine and dopamine receptor antagonist binding in the rat ventral and dorsal striatum: lack of changes after a neonatal bilateral lesion of the

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481