


SCH-23985

Chemical Properties

CAS No. : 129666-34-8
Formula: C₁₈H₂₀ClNO
Molecular Weight: 301.81
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	SCH-23985 is an antagonist of D1 dopamine receptor.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3133 mL	16.5667 mL	33.1334 mL
5 mM	0.6627 mL	3.3133 mL	6.6267 mL
10 mM	0.3313 mL	1.6567 mL	3.3133 mL
50 mM	0.0663 mL	0.3313 mL	0.6627 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481