


SERT inhibitor 69419

Chemical Properties

CAS No. : 305374-42-9
Formula: C₂₅H₁₆N₂O₆S
Molecular Weight: 472.47
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	SERT inhibitor 69419 inhibits serotonin reuptake.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1165 mL	10.5827 mL	21.1654 mL
5 mM	0.4233 mL	2.1165 mL	4.2331 mL
10 mM	0.2117 mL	1.0583 mL	2.1165 mL
50 mM	0.0423 mL	0.2117 mL	0.4233 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481