


Tigolaner

Chemical Properties

CAS No. : 1621436-41-6
Formula: C₂₁H₁₃ClF₈N₆O
Molecular Weight: 552.81
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Tigolaner is a antiparasitic drug candidate for veterinary use.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8089 mL	9.0447 mL	18.0894 mL
5 mM	0.3618 mL	1.8089 mL	3.6179 mL
10 mM	0.1809 mL	0.9045 mL	1.8089 mL
50 mM	0.0362 mL	0.1809 mL	0.3618 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481