# Data Sheet (Cat.No.T28997)



#### TP-271

### **Chemical Properties**

CAS No.: 1207284-17-0

Formula: C27H31FN4O8

Molecular Weight: 558.56

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

## **Biological Description**

Description TP-271, a tetracyclines antibiotic, is used potentially for the treatment of community-

acquired bacterial pneumonia (CABP).

### **Preparing Stock Solutions**

	1mg	5mg	10mg
1 mM	1.7903 mL	8.9516 mL	17.9032 mL
5 mM	0.3581 mL	1.7903 mL	3.5806 mL
10 mM	0.179 mL	0.8952 m <b></b> €	1.7903 mL
50 mM	0.0358 mL	0.179 mL	0.3581 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com