


3-Diazirinophenylretinal

Chemical Properties

CAS No. : 83374-72-5
Formula: C₁₈H₁₈N₂O
Molecular Weight: 278.35
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	3-Diazirinophenylretinal is a bioactive chemical.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5926 mL	17.963 mL	35.926 mL
5 mM	0.7185 mL	3.5926 mL	7.1852 mL
10 mM	0.3593 mL	1.7963 mL	3.5926 mL
50 mM	0.0719 mL	0.3593 mL	0.7185 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481