


4-Hydroxytoremifene

Chemical Properties

CAS No. : 110503-62-3
Formula: C₂₆H₂₈ClNO₂
Molecular Weight: 421.96
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	4-Hydroxytoremifene is a selective estrogen receptor modulator toremifene active metabolite.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3699 mL	11.8495 mL	23.6989 mL
5 mM	0.474 mL	2.3699 mL	4.7398 mL
10 mM	0.237 mL	1.1849 mL	2.3699 mL
50 mM	0.0474 mL	0.237 mL	0.474 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Edavana VK, Dhakal IB, Yu X, Williams S, Kadlubar S. Sulfation of 4-hydroxy toremifene: individual variability, isoform specificity, and contribution to toremifene pharmacogenomics. Drug Metab Dispos. 2012 Jun;40(6):1210-

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481