


A61603

Chemical Properties

CAS No. : 750531-54-5
Formula: C₁₄H₁₉N₃O₃S
Molecular Weight: 309.38
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	A61603 is an alpha-1A adrenergic receptor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2323 mL	16.1614 mL	32.3227 mL
5 mM	0.6465 mL	3.2323 mL	6.4645 mL
10 mM	0.3232 mL	1.6161 mL	3.2323 mL
50 mM	0.0646 mL	0.3232 mL	0.6465 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Skelly MJ, Ariwodola OJ, Weiner JL. Fear conditioning selectively disrupts noradrenergic facilitation of GABAergic inhibition in the basolateral amygdala. *Neuropharmacology*. 2017 Feb;113(Pt A):231-240. doi: 10.1016/j.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481