


AI 3-63044

Chemical Properties

CAS No. : 1092-03-1
Formula: C₁₈H₂₃Cl₂N₃
Molecular Weight: 352.3
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	AI 3-63044 is a bioactive chemical.
-------------	-------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8385 mL	14.1924 mL	28.3849 mL
5 mM	0.5677 mL	2.8385 mL	5.677 mL
10 mM	0.2838 mL	1.4192 mL	2.8385 mL
50 mM	0.0568 mL	0.2838 mL	0.5677 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481