


3-Nitropropanoic acid

Chemical Properties

CAS No. : 504-88-1
Formula: C₃H₅NO₄
Molecular Weight: 119.08
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	3-Nitropropanoic acid (Bovinoicidin) is an irreversible inhibitor of succinate dehydrogenase.
Targets(IC50)	Apoptosis,Dehydrogenase,Antibacterial

Solubility Information

Solubility	DMSO: 125 mg/mL (1049.71 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.3977 mL	41.9886 mL	83.9772 mL
5 mM	1.6795 mL	8.3977 mL	16.7954 mL
10 mM	0.8398 mL	4.1989 mL	8.3977 mL
50 mM	0.168 mL	0.8398 mL	1.6795 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Gupta S, Sharma B. Eur J Pharmacol. 2014 Jun 5;732:111-22.

Xiong J, Wang L, Feng Y, et al.Geldanamycin confers fungicidal properties to azole by triggering the activation of

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481