


AS100

Chemical Properties

CAS No. : 860033-28-9
Formula: C₂₃H₂₄Cl₃N₇O₄
Molecular Weight: 568.84
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	AS100 is a bioactive chemical.
-------------	--------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.758 mL	8.7898 mL	17.5796 mL
5 mM	0.3516 mL	1.758 mL	3.5159 mL
10 mM	0.1758 mL	0.879 mL	1.758 mL
50 mM	0.0352 mL	0.1758 mL	0.3516 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Jain A, Rao N, Sharifi M, Bhatt N, Patel P, Nirmal D, Persaud JW, Nair DR. Evaluation of the point of care Afinion AS100 analyser in a community setting. Ann Clin Biochem. 2016 Jul 15. pii: 0004563216661737. [Epub ahead of

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481