


Asterinin E

Chemical Properties

CAS No. : 172548-92-4
Formula: C₂₆H₃₅N₅O₉
Molecular Weight: 561.592
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Asterinin E is a pentapeptide extracted from the roots of Aster tataricus.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7807 mL	8.9033 mL	17.8066 mL
5 mM	0.3561 mL	1.7807 mL	3.5613 mL
10 mM	0.1781 mL	0.8903 mL	1.7807 mL
50 mM	0.0356 mL	0.1781 mL	0.3561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Cheng DL, Shao Y, Zhao K, Hartmann R, Roeder E. Pentapeptides from the roots of Aster tataricus. Pharmazie. 1996 Mar;51(3):185-6. PubMed PMID: 8900872.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481