


Bemarinone

Chemical Properties

CAS No. : 16535-98-1
Formula: C₁₂H₁₃NO₃
Molecular Weight: 219.24
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Bemarinone is a clinical cardiotonic agent.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5612 mL	22.8061 mL	45.6121 mL
5 mM	0.9122 mL	4.5612 mL	9.1224 mL
10 mM	0.4561 mL	2.2806 mL	4.5612 mL
50 mM	0.0912 mL	0.4561 mL	0.9122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Kanojia RM, Press JB, Lever OW Jr, Williams L, McNally JJ, Tobia AJ, Falotico R, Moore JB Jr. Cardiotoxic agents. Synthesis and inotropic activity of a series of isoquinolin-3-ol derivatives. J Med Chem. 1988 Jul;31(7):1363-8.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481