


Bolandiol dipropionate

Chemical Properties

CAS No. : 1986-53-4
Formula: C₂₄H₃₆O₄
Molecular Weight: 388.54
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Bolandiol dipropionate is an Anabolic.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5737 mL	12.8687 mL	25.7374 mL
5 mM	0.5147 mL	2.5737 mL	5.1475 mL
10 mM	0.2574 mL	1.2869 mL	2.5737 mL
50 mM	0.0515 mL	0.2574 mL	0.5147 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481