


C.I. 59830

Chemical Properties

CAS No. : 25704-81-8
Formula: C₃₆H₁₈Br₂O₄
Molecular Weight: 674.344
Appearance: no data available
Storage: keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	C.I. 59830 is a green dye.
-------------	----------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4829 mL	7.4147 mL	14.8293 mL
5 mM	0.2966 mL	1.4829 mL	2.9659 mL
10 mM	0.1483 mL	0.7415 mL	1.4829 mL
50 mM	0.0297 mL	0.1483 mL	0.2966 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481