

CN009543V

Chemical Properties


CAS No. : 375826-84-9

Formula: C₁₂H₁₂N₄O₆S

Molecular Weight: 340.31

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	CN009543V is an enhancer of tyrosine phosphorylation of EGFR via downstream signaling pathways.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9385 mL	14.6925 mL	29.385 mL
5 mM	0.5877 mL	2.9385 mL	5.877 mL
10 mM	0.2938 mL	1.4692 mL	2.9385 mL
50 mM	0.0588 mL	0.2938 mL	0.5877 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481