


Dihydroagathic acid

Chemical Properties

CAS No. : 5956-15-0
Formula: C₂₀H₃₂O₄
Molecular Weight: 336.47
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Dihydroagathic acid is a metabolite of isocupressic acid.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.972 mL	14.8602 mL	29.7203 mL
5 mM	0.5944 mL	2.972 mL	5.9441 mL
10 mM	0.2972 mL	1.486 mL	2.972 mL
50 mM	0.0594 mL	0.2972 mL	0.5944 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Snider DB, Gardner DR, Janke BH, Ensley SM. Pine needle abortion biomarker detected in bovine fetal fluids. J Vet Diagn Invest. 2015 Jan;27(1):74-9. doi: 10.1177/1040638714554443. Epub 2014 Nov 26. PubMed PMID: 25428187.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481