


FA-alkyne

Chemical Properties

CAS No. : 671188-20-8
Formula: C₇H₇NO₃
Molecular Weight: 153.137
Appearance: no data available
Storage: store at low temperature
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	FA-alkyne is a novel clickable fumarate-competitive chemoproteomic probe
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.530 mL	32.6499 mL	65.2997 mL
5 mM	1.306 mL	6.530 mL	13.0599 mL
10 mM	0.653 mL	3.265 mL	6.530 mL
50 mM	0.1306 mL	0.653 mL	1.306 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481