


## Fungiqua

### Chemical Properties

CAS No. : 12224-02-1  
Formula: C<sub>40</sub>H<sub>44</sub>N<sub>12</sub>NaO<sub>16</sub>S<sub>4</sub>  
Molecular Weight: 1100.09  
Appearance: no data available  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



### Biological Description

|             |                                   |
|-------------|-----------------------------------|
| Description | Fungiqua is a bioactive chemical. |
|-------------|-----------------------------------|

### Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg      |
|-------|-----------|-----------|-----------|
| 1 mM  | 0.909 mL  | 4.5451 mL | 9.0902 mL |
| 5 mM  | 0.1818 mL | 0.909 mL  | 1.818 mL  |
| 10 mM | 0.0909 mL | 0.4545 mL | 0.909 mL  |
| 50 mM | 0.0182 mL | 0.0909 mL | 0.1818 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

### Reference

Yamada S, Okudaira H, Hayashi E. An alteration in sensitivity to cholinergic agents on guinea-pig ilea and atria after repeated administration of an organophosphate and an antagonism by a carbamate. Arch Int Pharmacodyn

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481