


Gitoformate

Chemical Properties

CAS No. : 10176-39-3
Formula: C₄₆H₆₄O₁₉
Molecular Weight: 920.999
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Gitoformate is used to treat kidney and heart insufficiency,
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0858 mL	5.4289 mL	10.8578 mL
5 mM	0.2172 mL	1.0858 mL	2.1716 mL
10 mM	0.1086 mL	0.5429 mL	1.0858 mL
50 mM	0.0217 mL	0.1086 mL	0.2172 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Kubin S, Woodcock BG, Menke G, Rietbrock N. [Postprandial delay of drug absorption in a gitoformate model]. Arzneimittelforschung. 1987 Dec;37(12):1376-8. German. PubMed PMID: 3449065.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481