


GYKI-T412

Chemical Properties

CAS No. : 991-70-8
Formula: C₃₁H₃₇N₃O₅S
Molecular Weight: 563.71
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	GYKI-T412 is small molecule agent with potential bioactivity. Detail information has not been published.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.774 mL	8.8698 mL	17.7396 mL
5 mM	0.3548 mL	1.774 mL	3.5479 mL
10 mM	0.1774 mL	0.887 mL	1.774 mL
50 mM	0.0355 mL	0.1774 mL	0.3548 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481