


IMM-03

Chemical Properties

CAS No. : T32149
Formula: C₁₇H₂₇N₃S₂
Molecular Weight: 337.54
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	IMM-03 is a negative control for IMM-01 and IMM-02.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9626 mL	14.8131 mL	29.6261 mL
5 mM	0.5925 mL	2.9626 mL	5.9252 mL
10 mM	0.2963 mL	1.4813 mL	2.9626 mL
50 mM	0.0593 mL	0.2963 mL	0.5925 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481