


J 2644

Chemical Properties

CAS No. : 71712-03-3
Formula: C₂₂H₃₀O₂
Molecular Weight: 326.47
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	J 2644 is a bioactive chemical.
-------------	---------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0631 mL	15.3153 mL	30.6307 mL
5 mM	0.6126 mL	3.0631 mL	6.1261 mL
10 mM	0.3063 mL	1.5315 mL	3.0631 mL
50 mM	0.0613 mL	0.3063 mL	0.6126 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Tait A, Parenti C, Zanolini P, Veneri C, Truzzi C, Brandoli C, Baraldi M, di Bella M. Synthesis and antiinflammatory activity of 2,6-bis(1,1-dimethylethyl)phenol derivatives. *Farmaco*. 1993 Oct;48(10):1463-73. PubMed PMID:

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481