


K 13-004

Chemical Properties

CAS No. : 96609-38-0
Formula: C₂₂H₂₇NO₅
Molecular Weight: 385.45
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	K 13-004 is a lipid-lowering agent.
-------------	-------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5944 mL	12.9719 mL	25.9437 mL
5 mM	0.5189 mL	2.5944 mL	5.1887 mL
10 mM	0.2594 mL	1.2972 mL	2.5944 mL
50 mM	0.0519 mL	0.2594 mL	0.5189 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Schliack M, Löser R, Seibel K. Circadian variation of the hypocholesterolemic effect of K13.004 in rats. Atherosclerosis. 1987 Apr;64(2-3):163-6. PubMed PMID: 3606713.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481