


L 345

Chemical Properties

CAS No. : 121809-82-3
Formula: C₁₇H₁₅Cl₃N₂O₄
Molecular Weight: 417.67
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L 345 is a bioactive chemical.
-------------	--------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3942 mL	11.9712 mL	23.9423 mL
5 mM	0.4788 mL	2.3942 mL	4.7885 mL
10 mM	0.2394 mL	1.1971 mL	2.3942 mL
50 mM	0.0479 mL	0.2394 mL	0.4788 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Horjen AW, Seljeflot I, Berge T, Smith P, Arnesen H, Tveit A. Effect of sinus rhythm restoration on markers of thrombin generation in atrial fibrillation. Thromb J. 2017 Dec 28;15:30. doi: 10.1186/s12959-017-0153-1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481