


L 644257

Chemical Properties

CAS No. : 68354-99-4
Formula: C₃₉H₆₉NO₅
Molecular Weight: 664.03
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L 644257 is a synthetic glycolipid for study liposomes in vivo.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.506 mL	7.5298 mL	15.0596 mL
5 mM	0.3012 mL	1.506 mL	3.0119 mL
10 mM	0.1506 mL	0.753 mL	1.506 mL
50 mM	0.0301 mL	0.1506 mL	0.3012 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Ponpipom MM, Hagmann WK, O'Grady LA, Jackson JJ, Wood DD, Zweerink HJ. Glycolipids as host resistance stimulators. J Med Chem. 1990 Feb;33(2):861-7. PubMed PMID: 2299649.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481