


L 686741

Chemical Properties

CAS No. : 133768-71-5
Formula: C₃₁H₃₀Cl₂N₂O₄S
Molecular Weight: 597.55
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L 686741 is a bioactive chemical.
-------------	-----------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6735 mL	8.3675 mL	16.735 mL
5 mM	0.3347 mL	1.6735 mL	3.347 mL
10 mM	0.1674 mL	0.8368 mL	1.6735 mL
50 mM	0.0335 mL	0.1674 mL	0.3347 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Grossman SJ, DeLuca JG, Zamboni RJ, Keenan KP, Patrick DH, Herold EG, van Zwieten MJ, Zacchei AG.
Enantioselective induction of peroxisomal proliferation in CD-1 mice by leukotriene antagonists. Toxicol Appl

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481