


L023103

Chemical Properties

CAS No. : 1222539-85-6
Formula: C30H33ClN6O5S
Molecular Weight: 625.14
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L023103 is a bioactive chemical.
-------------	----------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5996 mL	7.9982 mL	15.9964 mL
5 mM	0.3199 mL	1.5996 mL	3.1993 mL
10 mM	0.160 mL	0.7998 mL	1.5996 mL
50 mM	0.032 mL	0.160 mL	0.3199 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481