


Laserolide

Chemical Properties

CAS No. : 26560-24-7
Formula: C₂₂H₃₀O₆
Molecular Weight: 390.476
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Laserolide is a sesquiterpene lactone.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.561 mL	12.8048 mL	25.6095 mL
5 mM	0.5122 mL	2.561 mL	5.1219 mL
10 mM	0.2561 mL	1.2805 mL	2.561 mL
50 mM	0.0512 mL	0.2561 mL	0.5122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Harmatha J, Buděšínský M, Vokáč K, Kostecká P, Kmoníčková E, Zídek Z. Trilobolide and related sesquiterpene lactones from *Laser trilobum* possessing immunobiological properties. *Fitoterapia*. 2013 Sep;89:157-66. doi:

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481