


L-Flamprop-isopropyl

Chemical Properties

CAS No. : 57973-67-8
Formula: C₁₉H₁₉ClFNO₃
Molecular Weight: 363.81
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L-Flamprop-isopropyl is a bioactive chemical.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7487 mL	13.7434 mL	27.4869 mL
5 mM	0.5497 mL	2.7487 mL	5.4974 mL
10 mM	0.2749 mL	1.3743 mL	2.7487 mL
50 mM	0.055 mL	0.2749 mL	0.5497 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Hakme E, Lozano A, Gómez-Ramos MM, Hernando MD, Fernández-Alba AR. Non-target evaluation of contaminants in honey bees and pollen samples by gas chromatography time-of-flight mass spectrometry.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481