


LY 117413

## Chemical Properties

CAS No. : 82970-70-5  
Formula: C<sub>22</sub>H<sub>29</sub>NO<sub>2</sub>  
Molecular Weight: 339.47  
Appearance: no data available  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



## Biological Description

|             |                                     |
|-------------|-------------------------------------|
| Description | LY 117413 is a bio-active chemical. |
|-------------|-------------------------------------|

## Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 2.9458 mL | 14.7288 mL | 29.4577 mL |
| 5 mM  | 0.5892 mL | 2.9458 mL  | 5.8915 mL  |
| 10 mM | 0.2946 mL | 1.4729 mL  | 2.9458 mL  |
| 50 mM | 0.0589 mL | 0.2946 mL  | 0.5892 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481