


## Methyl tetrahydropyranylraserpate

### Chemical Properties

CAS No. : 751-73-5  
Formula: C<sub>28</sub>H<sub>38</sub>N<sub>2</sub>O<sub>6</sub>  
Molecular Weight: 498.62  
Appearance: no data available  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



### Biological Description

|             |                                                     |
|-------------|-----------------------------------------------------|
| Description | Methyl tetrahydropyranylraserpate is a biochemical. |
|-------------|-----------------------------------------------------|

### Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 2.0055 mL | 10.0277 mL | 20.0554 mL |
| 5 mM  | 0.4011 mL | 2.0055 mL  | 4.0111 mL  |
| 10 mM | 0.2006 mL | 1.0028 mL  | 2.0055 mL  |
| 50 mM | 0.0401 mL | 0.2006 mL  | 0.4011 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

### Reference

Li X, Wang Z, Ning J, Gao M, Jiang W, Zhou Z, Li G. Preparation and characterization of a novel polyethyleneimine cation-modified persimmon tannin bioadsorbent for anionic dye adsorption. J Environ Manage. 2018 Mar 31;217:

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481