


MK-1-907

Chemical Properties

CAS No. : 78150-06-8
Formula: C₂₅H₃₂N₂O₂
Molecular Weight: 392.53
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	MK-1-907 is a potential bioactive agent. No further details yet.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5476 mL	12.7379 mL	25.4758 mL
5 mM	0.5095 mL	2.5476 mL	5.0952 mL
10 mM	0.2548 mL	1.2738 mL	2.5476 mL
50 mM	0.051 mL	0.2548 mL	0.5095 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481