


Nafetolol

Chemical Properties

CAS No. :	42050-23-7
Formula:	C ₁₉ H ₂₉ NO ₃
Molecular Weight:	319.44
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Nafetolol (K 5407) is a biochemical reagent that can be used to synthesize other compounds.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 45 mg/mL (140.87 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1305 mL	15.6524 mL	31.3048 mL
5 mM	0.6261 mL	3.1305 mL	6.261 mL
10 mM	0.313 mL	1.5652 mL	3.1305 mL
50 mM	0.0626 mL	0.313 mL	0.6261 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Goldaniga G, et al. Pharmacokinetics and metabolism of a new beta-adrenergic blocking agent, the 1, ter-butyl-amino-3-(1,2,3,4-tetrahydro-1,4-ethano-8-hydroxy-5-naphthoxy)-2-propanol (K 5407). Eur J Drug Metab

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481