


Neospiran

Chemical Properties

CAS No. : 83-81-8
Formula: C₁₆H₂₄N₂O₂
Molecular Weight: 276.37
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Neospiran can be used in intracardial injections.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6183 mL	18.0917 mL	36.1834 mL
5 mM	0.7237 mL	3.6183 mL	7.2367 mL
10 mM	0.3618 mL	1.8092 mL	3.6183 mL
50 mM	0.0724 mL	0.3618 mL	0.7237 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

PFEIFER W. [Value of intracardial neospiran injections]. Zentralbl Chir. 1952;77(18):773-5. Undetermined Language. PubMed PMID: 14959358.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481