


PF-06260414

Chemical Properties

CAS No. : 1612755-71-1
Formula: C₁₄H₁₄N₄O₂S
Molecular Weight: 302.35
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PF-06260414 is a selective androgen receptor modulator or SARM used to treat muscle weakness.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3074 mL	16.5371 mL	33.0743 mL
5 mM	0.6615 mL	3.3074 mL	6.6149 mL
10 mM	0.3307 mL	1.6537 mL	3.3074 mL
50 mM	0.0661 mL	0.3307 mL	0.6615 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481