


Phenarsone sulfoxylate

Chemical Properties

CAS No. : 535-51-3
Formula: C₇H₁₀AsNNa₂O₆S
Molecular Weight: 357.12
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Phenarsone sulfoxylate is a biochemical.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8002 mL	14.0009 mL	28.0018 mL
5 mM	0.560 mL	2.8002 mL	5.6004 mL
10 mM	0.280 mL	1.4001 mL	2.8002 mL
50 mM	0.056 mL	0.280 mL	0.560 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481