


Phenol, m-(1,4-dimethyl-3-propyl-3-pyrrolidinyl)-

Chemical Properties

CAS No. : 1505-23-3
Formula: C₁₅H₂₃NO
Molecular Weight: 233.35
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Phenol, m-(1,4-dimethyl-3-propyl-3-pyrrolidinyl)- is a bioactive chemical.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2854 mL	21.427 mL	42.8541 mL
5 mM	0.8571 mL	4.2854 mL	8.5708 mL
10 mM	0.4285 mL	2.1427 mL	4.2854 mL
50 mM	0.0857 mL	0.4285 mL	0.8571 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481