


Pteroyltyrosine

Chemical Properties

CAS No. : 76910-23-1
Formula: C₂₃H₂₁N₇O₅
Molecular Weight: 475.46
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Pteroyltyrosine is an analog of tyrosyl folate used in folate radioassay.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1032 mL	10.5161 mL	21.0323 mL
5 mM	0.4206 mL	2.1032 mL	4.2065 mL
10 mM	0.2103 mL	1.0516 mL	2.1032 mL
50 mM	0.0421 mL	0.2103 mL	0.4206 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Farina PR, Grattan JA. Radioiodinated pteroyltyrosine: a novel analog for folate radioassay. Anal Biochem. 1981 May 1;113(1):124-9. PubMed PMID: 7270877.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481