


Retroprogesterone

Chemical Properties

CAS No. :	2755-10-4
Formula:	C ₂₁ H ₃₀ O ₂
Molecular Weight:	314.46
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Retroprogesterone, a progestin, may be used for menstrual irregularities in adolescence, for the maintenance of threatened pregnancy, and possibly as an ovulation stimulant.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1801 mL	15.9003 mL	31.8005 mL
5 mM	0.636 mL	3.1801 mL	6.3601 mL
10 mM	0.318 mL	1.590 mL	3.1801 mL
50 mM	0.0636 mL	0.318 mL	0.636 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

AYDAR CK, GREENBLATT RB. 6-DEHYDRO-RETROPROGESTERONE (DUPHASTON) AN INTERESTING PROGESTERONE-LIKE COMPOUND. Int J Fertil. 1964 Oct-Dec;9:585-95. PubMed PMID: 14215356.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481