


Riseptin

Chemical Properties

CAS No. : 102-30-7
Formula: C₂₁H₃₆Cl₃N
Molecular Weight: 408.88
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Riseptin is a disinfectant.
-------------	-----------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4457 mL	12.2285 mL	24.4571 mL
5 mM	0.4891 mL	2.4457 mL	4.8914 mL
10 mM	0.2446 mL	1.2229 mL	2.4457 mL
50 mM	0.0489 mL	0.2446 mL	0.4891 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

IMHOLZ G. [Clinical usefulness of riseptin for disinfection of the hands]. Med Klin. 1957 Feb 8;52(6):224-5. German. PubMed PMID: 13418238.
KEUDELM J. [One year of use of riseptin in surgical hand disinfection]. Dtsch Med

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481