


Sarisan

Chemical Properties

CAS No. : 18607-93-7
 Formula: C₁₁H₁₂O₃
 Molecular Weight: 192.2112
 Appearance: no data available
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Sarisan is a biochemical.
-------------	---------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2026 mL	26.0132 mL	52.0264 mL
5 mM	1.0405 mL	5.2026 mL	10.4053 mL
10 mM	0.5203 mL	2.6013 mL	5.2026 mL
50 mM	0.1041 mL	0.5203 mL	1.0405 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Zhang ZW, Liu GX, Xie DM, Tian F, Jia YK, Xu F, Shang MY, Wang X, Cai SQ. [Quantitative determination of seven major absorbed volatile constituents in mice brain, liver and blood after intragastric administration of Asari Radix

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481