


SR 42128

Chemical Properties

CAS No. : 98092-14-9
Formula: C₃₉H₆₅N₅O₉
Molecular Weight: 747.975
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	SR 42128 is arginine salt.
-------------	----------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3369 mL	6.6847 mL	13.3693 mL
5 mM	0.2674 mL	1.3369 mL	2.6739 mL
10 mM	0.1337 mL	0.6685 mL	1.3369 mL
50 mM	0.0267 mL	0.1337 mL	0.2674 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Delie F, Couvreur P, Nisato D, Michel JB, Puisieux F, Letourneux Y. Synthesis and in vitro study of a diglyceride prodrug of a peptide. Pharm Res. 1994 Aug;11(8):1082-7. PubMed PMID: 7971705.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481