


Synstab A

Chemical Properties

CAS No. : 257612-36-5
Formula: C₁₅H₁₃BrCl₃N₃O₃S
Molecular Weight: 501.61
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Synstab A is a microtubule stabilizer.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9936 mL	9.9679 mL	19.9358 mL
5 mM	0.3987 mL	1.9936 mL	3.9872 mL
10 mM	0.1994 mL	0.9968 mL	1.9936 mL
50 mM	0.0399 mL	0.1994 mL	0.3987 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481