


Taltsv

Chemical Properties

CAS No. : 71487-43-9
Formula: C₃₁H₄₉N₇O₁₁
Molecular Weight: 695.76
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Taltsv is a bioactive peptide.
-------------	--------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4373 mL	7.1864 mL	14.3728 mL
5 mM	0.2875 mL	1.4373 mL	2.8746 mL
10 mM	0.1437 mL	0.7186 mL	1.4373 mL
50 mM	0.0287 mL	0.1437 mL	0.2875 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Palamarczyk G, Lehle L, Mankowski T, Chojnacki T, Tanner W. Specificity of solubilized yeast glycosyl transferases for polyprenyl derivatives. Eur J Biochem. 1980 Apr;105(3):517-23. PubMed PMID: 6445267.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481