


Telmapitant

Chemical Properties

CAS No. : 552292-58-7
Formula: C₂₄H₂₃F₆N₃O₃
Molecular Weight: 515.45
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Telmapitant is an antagonist of the neurokinin NK1 receptor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9401 mL	9.7003 mL	19.4005 mL
5 mM	0.388 mL	1.9401 mL	3.8801 mL
10 mM	0.194 mL	0.970 mL	1.9401 mL
50 mM	0.0388 mL	0.194 mL	0.388 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481